

ELISA Assay Kit for Anti-PD-1 h-mAb in Mouse Serum

Pack Size: 96 tests/480 tests

Catalog Number: EPM-V1

IMPORTANT: Please carefully read this manual before performing your experiment.

For Research Use Only. Not For Use In Diagnostic or Therapeutic Procedures

www.acrobiosystems.com



Intended Use

The enclosed kit is intended for research use only. The products are not intended for treatment and/or diagnosis of diseases.

Background

How a therapeutic antibody is metabolized in the body is pertinently relevant to its efficacy. Therefore, pharmacokinetics study is an important part of the drug development.

PD-1 is probably the most targeted molecule in today's pharmaceutical industry, thanking to the clinical success of Keytruda and Opdivo. Many investigational anti-PD-1 mAbs are being developed. There is a growing need for a standard assay that can be used to facilitate the study of their pharmacokinetics.

Assay Principles

The enclosed ELISA assay kit for anti-PD-1 h-mAb in mouse serum is based on an enzyme immunoassay (ELISA) between recombinant PD-1 protein (ECD) and the biotinylated anti-PD-1 monoclonal antibody. The method employs the principle of competitive ELISA, allowing quantification of the monoclonal anti-PD-1 antibodies in mouse serum.

The assay only involves the following steps:

a) Coat the plate with human PD-1;

b) Add the mixture of your sample and the biotinylated anti-PD-1 antibody;

c) Add Streptavidin-HRP followed by TMB or other colorimetric HRP substrate;

d) Record the OD readings and analyze the serum concentration.

Advantages

• Wide coverage

This kit has been tested for multiple investigational PD-1 antibodies, as well as the commercially available Keytruda and Opdivo. There is no requirement for the species and subtype of the antibody.

• Easy process

There is no longer a need to generate anti-drug antibodies for PK studies.

Clean result

The use of the biotinylated antibody in the competitive ELISA significantly alleviates the issue of background noise issue that often occurs in traditional ELISA method.

Materials Provided

Catalog	Components	Size (96 tests)	Size (480 tests)	Format	Storage
PK201-213	Human PD-1	15 µg	60 µg	Powder	- 20 ℃
PK202-213	Biotinylated anti- PD-1 antibody	10 µg	10 µg	Powder	-20 ℃
PK203-213	Streptavidin-HRP	5 µg	25 μg	Powder	-20℃, avoid light
PK204-213	Anti-PD-1 antibody	30 µg	30 µg	Powder	- 20 ℃

Tab.1 Materials Provided

 US and Canada:
 Tel: +1 800-810-0813

 Asia and Pacific:
 Tel: +86 400-685-2052

 Europe:
 Tel: +44 208-123-6158

Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110 Web: http://www.acrobiosystems.com

E-mail: order@acrobiosystems.com



Reconstitution

Reconstitute the provided lyophilized materials to stock solutions with PBS as recommended in **Tab.2.1 and Tab 2.2**, solubilize for 15 to 30 minutes at room temperature with occasional gentle mixing. **Avoid vigorous shaking or vortexing**.

The reconstituted stock solutions should be stored at -80 $^{\circ}$ C. **Avoid freeze-thaw cycles.** *Note: Streptavidin-HRP stock solution should be protected from light.*

Catalog	Components	Size	Stock Solution Con.	Reconstitution Buffer and Vol.
PK201-213	Human PD-1	15 µg	250 μg/mL	60 μL PBS
PK202-213	Biotinylated anti- PD-1 antibody	10 µg	50 μg/mL	200 μL PBS
PK203-213	Streptavidin-HRP	5 µg	50 μg/mL	100 μL PBS
PK204-213	Anti-PD-1 antibody	30 µg	200 μg/mL	150 μL PBS

Tab. 2.1 Reconstitution Method for 96 Tests

Tab. 2.2 Reconstitution Method for 480 Tests

Catalog	Components	Size	Stock Solution Con.	Reconstitution Buffer and Vol.
PK201-213	Human PD-1	60 µg	250 μg/mL	240 μL PBS
PK202-213	Biotinylated anti- PD-1 antibody	10 µg	50 μg/mL	200 μL PBS
PK203-213	Streptavidin-HRP	25 μg	50 μg/mL	500 μL PBS
PK204-213	Anti-PD-1 antibody	30 µg	200 μg/mL	150 μL PBS

Minimal Volume for Aliquot(s)

To avoid surface adsorp on loss and inactivation, the reconstituted protein must NOT be aliquoted to less than 10 μ g per vial.

Shipping and Storage

All components are shipped in lyophilized state at room temperature. This product is stable after storage at:

- 1) Room temperature (RT) for 1 month in lyophilized state;
- 2) -20 $^{\circ}$ C for 1 year in lyophilized state;
- 3) -80 $^\circ\!\mathrm{C}$ for 4 months under sterile conditions after reconstitution.



Materials or Instruments Required But Not Supplied (for 96 Tests)

Coating Buffer PBS (Phosphate Buffered Saline), pH7.4, 12 mL is sufficient.

Wash Buffer PBS with 0.05% (v/v) Tween-20 (PBST), 500 mL is sufficient.

Blocking Buffer Wash Buffer with 2% (w/v) bovine serum albumin (Sigma-Aldrich, Catalog # A4737), 35 mL is sufficient.

Dilution Buffer Wash Buffer with 0.5% (w/v) bovine serum albumin (Sigma-Aldrich, Catalog # A4737), 50 mL is sufficient.

Substrate Stock Solution 10 mg/mL TMB (Sigma-Aldrich, Catalog # 860336) in dimethyl sulfoxide (Sigma-Aldrich, Catalog # D8418), 1 mL is sufficient. **Protect from light.**

Substrate Dilution Buffer 50 mM disodium hydrogen phosphate (Na₂HPO₄) and 25 mM citric acid, adjust pH to 5.5 with 1 M sodium hydroxide (NaOH), 25 mL is sufficient.

TMB Substrate Working Solution

Dilute 250 μ L substrate stock solution in 25 mL substrate dilution buffer and add 12 μ L 5% H₂O₂ (pipette 10 μ L 30% H₂O₂ into 50 μ L distilled water), mix well.

Notes:

1) The TMB Substrate Working Solution should be freshly prepared and used within 15 minutes.

2) If you choose to use other commercially available ready-to-use TMB substrate solutions, you should follow the manufacturer's instruction.

Stop Solution 1 M sulfuric acid (aqueous), 6 mL is sufficient.

Therapeutic Anti-PD-1 Antibody to be Analyzed prepare by customer.

Serum Samples for Analyzed prepare by customer.

Mouse Serum Pre-dosed therapeutic serum or pooled by the equal amount of male and female serum as dilution matrix.

High Binding Surface 96-well Microplate, Clear Flat Bottom (Corning, Catalog # 42592)

Microplate Sealing Film (Sigma-Aldrich, Catalog # 28416050)

Pipettes and Pipette Tips

UV/Vis Microplate Spectrophotometer (absorbance 450 nm, correction wavelength set to 600 nm)

Note: Pilot experiment may be required to estimate the range of antibody concentration in samples.

Recommended Protocol (for 96 Tests)

1. Preparation

- 1.1. Reconstitute and store all reagents as recommended.
- Use the Therapeutic Anti-PD-1 Antibody to make 5 × Standards(STDs) Stock Solution with mouse serum. The 5 × STDs Stock Solution should include 0 μg/mL, 100 μg/mL, 50 μg/mL, 25 μg/mL, 12.5 μg/mL, 6.25 μg/mL, 3.125 μg/mL, 1.5625 μg/mL, 0.78125 μg/mL, 0.390625 μg/mL and 0.1953125 μg/mL of STDs respectively and stored frozen in single-use aliquots at -80°C until use.
- 1.3. Use the Therapeutic Anti-PD-1 Antibody to make 5 × Quality Controls(QCs) Stock Solution with mouse serum. The 5 × QCs Stock Solution should cover the quantification range of the standard curve (see below). Make independent dilutions from the antibody stock solution (e.g. include 25 µg/mL, 20 µg/mL, 5 µg/mL, 2 µg/mL and 0.78525 µg/mL of QCs) and stored frozen in single-use aliquots at -80°C until use.



1.4. Use the provided **Anti-PD-1 Antibody** as a reference if you need; dilute the antibody with mouse serum as recommended in 1.2.

2. Coating

- 2.1. Dilute human **PD-1 Stock Solution** (250 μg/mL) to 1 μg/mL with **Coating Buffer** to make human **PD-1 Working Solution**.
- 2.2. Please leave a couple of wells uncoated for STD-1 (No-Coating) and Only 10% Serum (No-Coating), respectively (Fig.1).
- 2.3. Add 100 μ L of human **PD-1 Working Solution** (1 μ g/mL) to each well, seal the plate with microplate sealing film and incubate overnight (or 16 hours) at 4 °C.

3. Washing

Remove the remaining solution by aspiration, add 300 μ L of **Wash buffer** to each well, gently tap the plate for 1 minute, remove any remaining **Wash Buffer** by aspirating or decanting, invert the plate and blot it against paper towels. **Repeat the wash step above for three times.**

4. Blocking

Add 300 μL **Blocking Buffer** to each well at 37 $^\circ\!C$ for 1.5 hours.

5. Washing

Repeat step 3. At meantime, you can start to prepare your samples.

6. Add Samples

- 6.1. Dilute the **Biotinylated Anti-PD-1 Antibody Stock Solution** (50 μg/mL) to 0.02 μg/mL with **Dilution Buffer** to make **Biotinylated Anti-PD-1 Antibody Working Solution**.
- 6.2. Dilute 5 × STDs and 5 × QCs with Dilution Buffer, then mixed with same volume Biotinylated Anti-PD-1 Antibody Working Solution. (For example: 22 μ L 5 × STDs + 88 μ L Dilution Buffer +110 μ L Biotinylated Anti-PD-1 Antibody Working Solution).
- 6.3. Dilute patient Samples with same volume Biotinylated Anti-PD-1 Antibody Working Solution.
- 6.4. For STD-1 (No-Coating), dilute 5×STD-1 (100 μg/mL) with Dilution Buffer, then mixed with same volume Biotinylated Anti-PD-1 Antibody Working Solution (For example: 22 μL 5×STD-1 + 88 μL Dilution Buffer +110 μL Biotinylated Anti-PD-1 Antibody Working Solution). Add 100 μL mixer to the wells.
- 6.5. For Only 10% Serum (No-Coating) and Only 10% Serum (Coating), please add 100 μL Dilution Buffer with 10% Mouse Serum (Fig. 1).
- 6.6. For Only Dilution Buffer (Coating), please add 100 μL Dilution Buffer (Fig. 1).
- 6.7. For all other wells, add 100 μ L mixed **STDs, QCs and Samples** respectively, seal the plate with microplate sealing film and incubate at 37 $^{\circ}$ C for 1 hour.

7. Washing

Repeat step 3.

8. Streptavidin-HRP Conjugating

- 8.1. Dilute Streptavidin-HRP Stock Solution (50 μ g/mL) to 0.1 μ g/mL with Dilution Buffer to make Streptavidin-HRP Working Solution.
- 8.2. For all wells, add 100 μ L **Streptavidin-HRP Working Solution**, seal the plate with microplate sealing film and incubate at 37 $^{\circ}$ C for 1 hour, avoid light.

Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110



9. Washing

Repeat step 3.

10. TMB Substrate Reaction

Add 200 μ L TMB Substrate Working Solution to each well. Incubated 37 $^{\circ}$ C for 20 minutes, avoid light.

11. Termination

Add 50 μ L **Stop Solution** to each well, and tap the plate gently for 3 minutes to allow thorough mixing. *Note:* The color in the wells should change from blue to yellow.

12. Data Recording

Read the absorbance at 450 nm using UV/Vis microplate spectrophotometer.

Note: The plate may be read at 600 nm without adding 1M sulfuric acid, but the Signal-to-Background ratio may be reduced.



Fig.1 Plate Layout

 US and Canada:
 Tel: +1 800-810-0813

 Asia and Pacific:
 Tel: +86 400-685-2052

 Europe:
 Tel: +44 208-123-6158

Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110



Tab.3 Assay Protocol

Steps Code	Steps	Reagents & Instruments	R e a c ti o n Conditions	Samples/ STDs/QCs	STD-1 (No- Coating)	Only 10% Serum (No- Coating)	Only 10% Serum (Coating)	Only Dilution Buffer (Coating)
1	Preparation	N/A	N/A	N/A	N/A	N/A	N/A	N/A
2	Coating	Human PD-1 working solution	4℃ for overnight	100 μL	-	-	100 μL	100 µL
3	Washing	Wash Buffer	Wash for 3 times	300 μL	300 μL	300 μL	300 μL	300 μL
4	Blocking	Blocking Buffer	37℃ for 1.5 hours	300 μL	300 μL	300 μL	300 μL	300 μL
5	Washing	Wash Buffer	Wash for 3 times	-	100 μL	-	-	-
6	Add Samples	Mixed samples/ Mixed $1 \times STDs/$ Mixed $1 \times QCs$	-	100 µL	-	-	-	
		Mixed 1×STD-1	-	-	100 μL	-	-	-
		10% Serum	-	-	-	100 μL	100 μL	-
		Dilution Buffer	-	-	-	-	-	100 µL
7	Washing	Wash Buffer	Wash for 3 times	300 μL	300 μL	300 μL	300 μL	300 μL
8	Streptavidin -HRP Labeling	Streptavidin -HRP Working Solution	37℃ for 1 hour	100 µL	100 μL	100 µL	100 µL	100 µL
9	Washing	Wash Buffer	Wash for 3 times	300 μL	300 μL	300 μL	300 μL	300 μL
10	TMB Substrate Reaction	TMB Substrate Working Solution	37℃ for 20 minutes	200 μL	200 μL	200 μL	200 μL	200 μL
11	Termination	Stop Solution	Mix by gentle tapping for 3 minutes	50 μL	50 μL	50 μL	50 μL	50 μL
12	Data Recording	UV/Vis spectrophotometer	Measure absorbance at 450 nm					

Note: It is recommended that all samples, controls and standards should be done in duplicates.

Method Verification

• PD-1: Biotinylated Anti-PD-1 Antibody Binding in the Absence of Competitors

Immobilized human **PD-1 Protein** at 1 µg/mL (100 µL/well) can bind the **Biotinylated Anti-PD-1 Antibody** with a linear range of 0.3-5 ng/mL when detected by **Streptavidin-HRP**. Background was subtracted from data points before curve fitting.





Fig.2 Binding of the biotinylated anti-PD-1 antibody to immobilized human PD-1 in a functional ELISA assay.

Sample Con. (µg/ml)	Mean Abs (OD450)
0.1	2.698
0.05	2.665
0.025	2.621
0.0125	2.616
0.00625	2.152
0.003125	1.269
0.0015625	0.648
0.00078125	0.341
0.000390625	0.160
0.0001953125	0.088
Blank	0.070

• Inhibition of PD-1: Biotinylated Anti-PD-1 Antibody Binding by An Anti-PD-1 Antibody

Serial dilutions of an **Anti-PD-1 Antibody** (Catalog # PK204-213) (1:2 serial dilutions, from 20 μ g/mL to 0.039 μ g/mL) was added into **PD-1: Biotinylated Anti-PD-1 Antibody** binding reactions. The assay was performed according to the above described protocol. Background was subtracted from data points prior to log transformation and curve fitting.



Fig.3 Inhibition of PD-1: biotinylated anti-PD-1 antibody binding by an anti-PD-1 antibody.

Sensitivity: The lowest detectable level that can be distinguished from the zero standard is 0.15625 μ g/mL. **Detection Range:** 0.15625-5 μ g/ml.

Precision:

Intra-assay CV: <20% for both calibration and QC points range 5- 0.15625 μg/mL for test run. Relax to CV < 25% and RE < 25% at the LLOQ (lower limit of quantification) and ULOQ (upper limit of quantification); for QC samples at high/low end use the same < 25% rule.

Inter-assay CV: < 20% for both calibration and QC points range 5- 0.15625 μ g/mL for test run. Relax to CV < 25% and RE < 25% at the LLOQ (lower limit of quantification) and ULOQ (upper limit of quantification); for QC samples at high/low end use the same < 25% rule.

Recovery: Recovery rate was found to be between 80-120% with normal mouse serum samples with known concentrations.

Note: All final data should be within the detection range as described above.

US and Canada: Tel: +1 800-810-0813 Asia and Pacific: Tel: +86 400-685-2052 Europe: Tel: +44 208-123-6158 Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110 Web: http://www.acrobiosystems.com

E-mail: order@acrobiosystems.com



• Quantitative Analysis of An Anti-PD-1 Antibody in Serum Samples

Serial dilutions of the **Therapeutic Anti-PD-1 Antibody** (1:2 serial dilutions, from 20 μ g/mL to 0.039 μ g/mL) was added into **PD-1: Biotinylated Anti-PD-1 Antibody** binding reactions. The assay was performed according to the above described protocol. Background was subtracted from data points prior to log transformation and curve fitting.



Therapeutic PD-1 Therapeutic PD-1 **Therapeutic PD-1** Pembrolizumab Nivolumab Antibody H, Name Antibody J, Antibody X, Human IgG4 Human IgG4 Human IgG4 Detection 0.03125-20 0.03125-20 0.03125-20 0.03125-20 0.03125-20 Rage (µg/ml) Sensitivity 0.15625 0.15625 0.15625 0.15625 0.15625 $(\mu g/ml)$ 87-115 %Recovery 90-105 94-109 89-112 97-109

Fig.4 Quantitative analysis five different therapeutic anti-PD-1 antibody in mouse serum samples.

Troubleshooting Guide

In case of a failed experiment, please check the expiration dates and status of the individual reagents, and make sure that all reagents have been reconstituted and stored as recommended. In addition, please make sure that all equipments are functioning properly.

Below is a list of common problems, and some tips that may help solve the problems and improve your results. *If you have any questions, please contact our technical support team at:TechSupport@acrobiosystems.com*

Tab.4	Troub	leshooting	Guide
-------	-------	------------	-------

Problem	Possible Cause	Solutions	
		Be sure the blocking step is performed.	
High background	Insufficient washing or blocking	 Increase number of washes and the volume Wash Buffer used. Increase Tween-20 concentration to 0.1% in Wash Buffer. Make sure Streptavidin-HRP is diluted in Blocking Buffer. 	
	Sample solvent contains inhibiting factors	 Run a negative control assay with the solvent alone. Maintain DMSO level at <1%. Increase protein incubation time. 	

US and Canada: Tel: +1 800-810-0813 Asia and Pacific: Tel: +86 400-685-2052 Europe: Tel: +44 208-123-6158 Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110



Tab. 4 Troubleshooting G	Guide (Continued)
--------------------------	-------------------

Problem	Possible Cause	Solutions	
	Contamination	- Make sure buffers and samples are prepared, used and stored correctly.	
	The TMB Substrate Working Solution is not fresh	- TMB Substrate Working Solution must be used within 15 minutes after preparation.	
	Inconsistent pipetting or dilution methods	 Make sure pipettors are functioning properly and use a multichannel pipettor if possible. Use master mixes to minimize errors. Run duplicates for all tests. 	
Colorimetric signal is erratic	TMB Substrate Working Solution is not completely mixed with the reaction solution	 Make sure that TMB Substrate Working Solution is adequately mixed with the reaction solution. 	
	Bubbles in the wells	- Tap plate gently to disperse bubbles.	
	Signal is too low	 The concentration of the samples should be adjusted to achieve optimal reading. Increase colorimetric HRP substrate incubation time 	
	Human PD-1, human anti-PD-1 antibody, or Streptavidin - HRP may have lost activity	 Make sure your proteins are aliquoted into single-use aliquots. Increase the time of reaction or increase the protein concentrationmay help in case the protein activity is decreased over time. 	
	Errors in instrument settings	- Please check instrument setting.	
Signal of positivecontrol is weak or abnormal	Substrate Stock Solution is outdated; Incubation temperature is incorrect; Incubation time is not sufficient	 Make sure the Substrate Stock Solution is working. Use proper incubation time and temperature. 	
	Insufficient mix	- Make sure the sample mixed sufficient before add to the plate	
	Pipetting errors	- Make sure that the pipette is calibrated and	
	Incomplete removal of residual buffers during previous steps	- Wells should appear dry after aspiration.	
Inadequate colorDevelopment	Problems with conjugate or color reagents	- Color should appear immediately after the reagent is added. Make sure no contamination or residual buffers in the wells before you start the color development process.	

Fax: +1 888-377-6111 Fax: +86 10-5885-0566 Fax: +44 203-355-6110 Web: http://www.acrobiosystems.com

E-mail: order@acrobiosystems.com