

Monkeypox virus (strain Zaire-96-I-16) A35R Protein, His Tag

Catalog # A3R-M52H3



Synonym

A35R (MPXV)

Source

Monkeypox virus (strain Zaire-96-I-16) A35R , His Tag(A3R-M52H3) is expressed from human 293 cells (HEK293). It contains AA Arg 58 - Thr 181 (Accession # [Q8V4U4](#)).
Predicted N-terminus: Arg 58

Molecular Characterization

A35R(Arg 58 - Thr 181)
Q8V4U4

Poly-his

This protein carries a polyhistidine tag at the C-terminus.
The protein has a calculated MW of 15.5 kDa. The protein migrates as 15 kDa and 17 kDa when calibrated against [Star Ribbon Pre-stained Protein Marker](#) under reducing (R) condition (SDS-PAGE) due to glycosylation.

Endotoxin

Less than 1.0 EU per µg by the LAL method / rFC method.

Purity

>90% as determined by SDS-PAGE.

Formulation

Lyophilized from 0.22 µm filtered solution in PBS, pH7.4 with trehalose as protectant.
Contact us for customized product form or formulation.

Reconstitution


Please see Certificate of Analysis for specific instructions.
For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Please avoid repeated freeze-thaw cycles.
This product is stable after storage at:

- 20°C to -70°C for 12 months in lyophilized state;
- 70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE




Monkeypox virus (strain Zaire-96-I-16) A35R , His Tag on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

Discounts, Gifts,
and more!

Monkeypox virus (strain Zaire-96-I-16) A35R Protein, His Tag

Catalog # A3R-M52H3



Immobilized Monkeypox virus (strain Zaire-96-I-16) A35R , His Tag (Cat. No. A3R-M52H3) at 5 µg/mL (100 µL/well) can bind Anti-Monkeypox-A35R-Antibody with a linear range of 0.1-3 ng/mL (QC tested).

Background

Monkeypox is a rare zoonosis caused by monkeypox virus, which has become the most serious orthpoxvirus and consists of complex double stranded DNA. The cases are mostly in central and western Africa. The pathogenesis of monkeypox is that the virus invades the body from respiratory mucosa, multiplies in lymphocytes, and incurs into blood producing transient venereal toxemia. after the virus multiplies in cells, the cells can invade the blood and propagate to the skin of the whole body, causing lesions. The envelope glycoprotein A35R on the EV surface has been predicted to influence intercellular diffusion of virions.

Discounts, Gifts,
and more!

